


Multiple Gas Analyzer + Sulfur GC Systems



- **Detect Low Sulfur Concentrations**
- **Multiple Gas AND Sulfur Analysis in one unit**
- **Uses just 1 Gas Sampling Valve and 3 Columns**
- **Room temperature Silcosteel Sample Loop**
- **Built-in, “whisper quiet” Air Compressor**
- **4 channel PeakSimple Data System**
- **...on the compact 8610C chassis**

Many analysts require natural gas analysis for BTU value calculations or drilling and mudlogging applications. Frequently, sulfur compounds are also of interest.

Because low sulfur concentrations (<50ppm) are difficult to measure, SRI has enhanced our popular Multiple Gas Analyzer GCs to simultaneously monitor low levels of sulfur compounds. The additional hardware required is an FPD/FID detector, which selectively detects sulfur down to mid-ppb range, and a room temperature Silcosteel sample loop.

One reason sulfur is so difficult to measure is that it disappears on contact with hot stainless steel surfaces; even limited contact with a hot stainless steel sample loop will destroy any sulfur in the gas sample. Our solution is to use a Silcosteel-lined transfer line leading to a splitter, and a 60-meter thick film capillary column. While Teflon columns are sometimes also used for sulfur analysis, the natural gas analysis (using MoleSieve and SilicaGel) requires column temperatures of 250°C or higher. Since the sulfur column is located in the same column oven, it is essential to use a column like the 60-meter capillary which can tolerate the higher temperatures.

8610-0073 Multiple Gas Analyzer #1 + Sulfur GC with TCD, FID, and FPD/FID detectors, methanizer, built-in air compressor, 3 columns, and Silcosteel sample loop

8610-0273 Multiple Gas Analyzer #2 + Sulfur GC with TCD, FID-methanizer, and FPD/FID detectors, built-in air compressor, 3 columns, and Silcosteel sample loop



OPTIONS & UPGRADES: 6 channel USB PeakSimple Data System, Split/Splitless and PTV injectors, additional column(s), Gas Sampling Valve, Thermal Desorber
(VOLTAGE: for 110VAC, use “part number-1” [ex: 8610-0073-1] for 220VAC, use “part number-2”)